## **RIDASCREEN® EBV VCA**

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Enzyme immunoassay for the detection of  
antibodies (IgG, IgM) against capsid antigen of EBV

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RIDA® and RIDASCREEN®  
are registered trademarks of R-Biopharm GmbH  
Manufacturer: R-Biopharm GmbH, Darmstadt, Germany

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Contents

	page
1. Intended use.....	15
2. General.....	15
3. Test principle.....	16
4. Reagents provided.....	17
5. Materials required but not provided.....	18
6. Warnings and precautions for the users.....	19
7. Storage instructions.....	19
8. Indication of instability or deterioration of reagents.....	20
9. Specimen collection and storage.....	20
10. Test procedure.....	20
11. Analysis.....	23
12. Remarks about the test procedure and interpretation.....	25

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## RIDASCREEN® EBV VCA

Enzyme immunoassay for the detection of antibodies (IgG, IgM) against capsid antigen of EBV

### 1. Intended use

The RIDASCREEN® EBV VCA test is an enzyme immunoassay (EIA) for the detection of antibodies (IgG and IgM) against capsid antigen of Epstein-Barr virus (EBV) in human serum.

### 2. General

The Epstein-Barr virus (EBV) belongs to the family of  $\gamma$ -Herpes viruses. Its systematic name is Human Herpes virus 4. As with other Herpes viruses, the genome consists of double-stranded DNA. This is contained in an enclosed, icosahedric capsid. The target cells of the EBV are human B- and T-lymphocytes within which the virus initially reproduces itself lytically releasing new infectious virus particles at the same time. After this, the virus transforms into the latent stage typical for Herpes viruses. The virus remains in this state in some of the B cells. As well as changing into the latent stage, a small amount of the virus is excreted via the mucous membrane of the mouth throughout the life of the host.

The incubation time for the virus is approximately four to six weeks. In the case of children, the course of the infection often runs without any symptoms. For adolescents and adults, the virus develops the clinical symptoms of an infectious mononucleosis after the incubation period. This self-limiting disease, often lasting for weeks, is also known as Pfeiffer's disease or glandular fever. The symptoms are sore throat, fever and swollen lymph nodes which are often accompanied by watering eyes. After two or three weeks, another symptom which often appears is a swollen spleen and, in a quarter of patients, skin rashes. In rare cases, chronic permanent infection can develop and, while symptoms persist, the virus transforms into the latent stage to a small degree only or not at all. Furthermore,

multiplication occurs lytically. During this process the symptoms which appear are tiredness, exhaustion and swelling of the lymph nodes, but this is not as pronounced as it is during the acute stage of mononucleosis. In the case of patients with immune suppression, there is also the danger of a B-cell lymphoma. For organ transplant patients, the risk of this happening increases by up to 5% while in the case of AIDS patients this risk is as high as 60%.

In Eastern and South East Asia, EBV is also an important factor in the occurrence of nasopharyngeal carcinoma. At the rate of 20 cases per 100,000 each year, this is the most frequently occurring cancer in this region.

For the diagnosis, it is possible to distinguish between the ways in which the antibodies act against the different virus proteins: proteins from virus replication (EA = early antigen), structure proteins (VCA = virus capsid antigen) and nuclear antigens (EBNA = Epstein-Barr nuclear antigens). IgM, IgA and IgG antibodies acting against VCA are generally found as the symptoms start to occur. After a few months, IgM and IgA antibodies are no longer found and therefore serve as a marker for primary infection. IgG antibodies acting against the structure proteins remain for the rest of the patient's life. IgM and IgA antibodies acting against early antigens can also be detected early on and can also be used as a marker for primary infection. The formation of IgG antibodies acting against EA has a somewhat delayed start. Furthermore, in contrast to IgG antibodies acting against VCA, IgG antibodies acting against EA do not remain for the whole of the patient's life but deteriorate after a few months. After three to six months, with the transition of the virus into its latent stage, the formation of IgG antibodies against nuclear antigens starts to take place. Of importance are especially IgG antibodies against the EBNA1 protein. This anti-EBNA1 IgG serves as a marker for the end of the acute infection phase and remains in existence for the rest of the patient's life. By non-appearance the latent stage (chronic persistent infection), these antibodies are missing. In the case of patients with immune suppression, they may also decrease again and disappear altogether.

IgA antibodies are important for the diagnosis of nasopharynx carcinoma, which is rare in Europe. An increased titer against these antibodies after the acute infection phase is over may indicate the onset of a carcinoma. Timely diagnosis can then lead to its effective treatment.

### 3. Test principle

On the surface of the microtiter wells, purified capsid antigen of EBV are bound. Diluted serum samples and Controls are pipetted into the wells and incubated at 37 °C in a humid chamber. Present antibodies bind to the immobilized antigen. Unbound material is removed in a washing step.

In a second step, a peroxidase-conjugated anti-human antibody (anti-IgG or IgM) is added. After incubation, unbound conjugate is removed by washing. Substrate (H<sub>2</sub>O<sub>2</sub>/TMB) is added to the wells and incubated at 37 °C in a humid chamber. The enzyme bound in the wells converts the colorless Substrate to a blue color. Addition of Stop Solution converts the color from blue to yellow. The absorption is measured at 450 nm wavelength (reference wavelength ≥ 620 nm).

### 4. Reagents provided

#### 4.1. RIDASCREEN® EBV VCA IgG (Art. No. K 6721)

The reagents in one package are sufficient for 96 determinations. Each test kit contains:

- 1 x 12 Microwell Strips (breakable) with 8 wells each in a frame; coated with capsid antigen of EBV; in a resealable foil bag
- 1 x RIDASCREEN® SeroPP (110 ml); Sample Diluent, ready to use; contains 0.01% Thimerosal
- 1 x RIDASCREEN® SeroWP (100 ml, 10x conc.); Washing Buffer, contains 0.2% Bronidox-L
- 1 x Standard Control (2.5 ml); diluted human serum, ready to use; contains 0.01% Thimerosal
- 1 x Negative Control (1.2 ml); diluted human serum, ready to use; contains 0.01% Thimerosal
- 1 x RIDASCREEN® SeroG LD (12 ml); Anti-*human-IgG*-Conjugate, ready to use; peroxidase-conjugated antibody (goat)
- 1 x RIDASCREEN® SeroSC (12 ml); H<sub>2</sub>O<sub>2</sub>/tetramethylbenzidine (TMB), ready to use
- 1 x RIDASCREEN® SeroStopp (12 ml); 1 N sulfuric acid
- 1 x Instructions for use

## 4.2. RIDASCREEN® EBV VCA IgM (Art. No. K 6731)

The reagents in one package are sufficient for 96 determinations.  
Each test kit contains:

- 1 x 12 Microwell Strips (breakable) with 8 wells each in a frame;  
coated with capsid antigen of EBV;  
in a resealable foil bag
- 1 x RIDASCREEN® SeroPP (110 ml);  
Sample Diluent, ready to use;  
contains 0.01% Thimerosal
- 1 x RIDASCREEN® SeroWP (100 ml, 10x conc.);  
Washing Buffer, contains 0.2 % Bronidox-L
- 1 x Standard Control (2.5 ml);  
diluted human serum, ready to use;  
contains 0.01% Thimerosal
- 1 x Negative Control (1.2 ml);  
diluted human serum, ready to use;  
contains 0.01% Thimerosal
- 1 x RIDASCREEN® SeroM LD (12 ml);  
Anti-human-IgM-Conjugate, ready to use;  
peroxidase-conjugated antibody (goat)
- 1 x RIDASCREEN® SeroSC (12 ml);  
H<sub>2</sub>O<sub>2</sub>/tetramethylbenzidine (TMB), ready to use
- 1 x RIDASCREEN® SeroStopp (12 ml);  
1 N sulfuric acid
- 1 x Instructions for use

## 5. Reagents required but not provided

### 5.1 Reagents

- Distilled or deionized water

### 5.2 Accessories

- Moist chamber at 37 °C
- Test tubes
- Vortex mixer
- Micropipets for volumes of 10 - 100 µl and 100 - 1000 µl
- Measuring cylinder (1000 ml)
- Microplate washer or multichannel pipet
- Microplate reader (450 nm, optional reference wavelength ≥ 620 nm)
- Absorbent paper

## 6. Warnings and precautions for the users

The control sera (Standard Control and Negative Control) have been tested for HIV- and HCV-Ab as well as for HBsAg and were found to be negative. However, they as well as the patient samples should be considered potentially contagious and be treated with the necessary safety precautions.

The Standard Control and the Negative Control as well as the Sample Diluent contain 0.01 % Thimerosal. Contact with skin or mucous membranes must be avoided.

The Washing Buffer contains Bronidox-L. Contact with skin or mucous membranes must be avoided.

Hydrogen peroxide can cause cauterization. Handle with care!

The Stop Solution contains 1 N sulfuric acid. Avoid contact with skin and clothing!

All reagents and materials coming in contact with potential infectious specimens must be treated with disinfectants or autoclaved at 121 °C for at least one hour.

An exchange of the Controls between kits of different lot numbers is not possible.

## 7. Storage instructions

All reagents have to be stored at 2 - 8 °C and can be used up to the expiry date printed on the labels. Microbial contamination has to be avoided. A quality warranty cannot be given beyond the kit expiration date.

The diluted Washing Buffer (RIDASCREEN® SeroWP) can be used up to the expiry date printed on the labels if stored at 2 - 8 °C.

Allow reagents and Microwell Strips to get room temperature before use. To avoid moisture within the strips, do not take the strips out of the foil bag before having reached room temperature. The foil bag should be opened with a pair of scissors without detaching the fastener. Return any unused strips to the foil bag, reseal and store them directly at 2 - 8 °C.

It is absolutely necessary to avoid contamination of the Substrate (RIDASCREEN® SeroSC) with the conjugate solution, since this results in a coloration of the Substrate. In the same way, the colorless Substrate must be protected from exposure to direct light to avoid deterioration or coloration by autoxidation. If the Substrate has turned blue, the reagent should be discarded.

## 8. Indication of instability or deterioration of reagents

The following criteria may indicate a reagent deterioration:

- a turbidity or a blue coloration of the colorless Substrate prior to its use
- an absorbance value of the Negative Control at 450/620 nm > 0.3
- an absorbance value of the Positive Control at 450/620 nm out of the values given in the attached tables of values

## 9. Specimen collection and storage

The RIDASCREEN® EBV VCA EIA has been evaluated for the investigation of human serum samples. Repeated freezing and thawing of the samples as well as microbial contamination must be avoided. The application of heat treated, lipemic, hemolytic, icteric or turbid samples can lead to wrong results.

The sample material can be stored for up to 1 week at 2 - 8 °C if the test cannot be carried out immediately. A prolonged storage of samples is possible at -20 °C.

Diluted samples can be used for up to 7 hours if they have been stored at 2 - 8 °C.

## 10. Test procedure

### 10.1. Preliminary comments

The RIDASCREEN® EIAs for serological antibody detection contain some test independent reagents. These reagents are labeled RIDASCREEN® Sero (SeroWP, SeroPP, SeroG LD, SeroM LD, SeroSC, SeroStopp). They are not kit specific and can be changed between single test kits with correspondingly labeled reagents. But, this is not valid for the Controls.

Bring all reagents and the Microwell Strips to room temperature before use. Mix the reagents well before use. Reproducibility in any EIA depends on exact pipetting, the observance of incubation times and temperature and the consistency of wash sequences. A deviation from the given incubation times and temperatures leads to an excessive difference of the standard control serum to the given target value, which is no longer covered by the determined area of values. During the washing steps, take care that all wells are filled with buffer and that the liquid is completely removed from the wells. Do not allow Microwells to dry between steps.

Avoid direct sunlight during all incubations. Covering the microtiter plate is recommended.

Except the Washing Buffer, all reagents are ready to use.

### 10.2. Preparation of the Washing Buffer

1 part of the concentrated Washing Buffer (RIDASCREEN® SeroWP) is diluted with 9 parts of distilled water. Crystals in the buffer concentrate can be dissolved in a waterbath at 37 °C. Add 100 ml of the concentrated Washing Buffer to a 1000 ml graduated cylinder. Bring the final volume to 1000 ml with distilled or deionized water. The diluted Washing Buffer can be used up to the expiry date printed on the labels if stored at 2 - 8 °C.

### 10.3. Preparation of the samples

Before starting the test, serum samples have to be diluted 1:100 with the Sample Diluent (RIDASCREEN® SeroPP).

Dilution of serum samples (IgG and IgM determination):

e.g. 10 µl serum + 990 µl Sample Diluent

For IgM determination, sera should be treated with an IgG-adsorbent (e.g. RIDA® G-Adsorbens, Art. No. Z 0101/Z 0102) prior to the test, following that it should be brought to its final dilution of 1:100 with the Sample Diluent.

#### **Attention!**

**Negative Control and Standard Control are ready to use and must not be diluted or adsorbed.**

### 10.4. First incubation

After insertion of a sufficient number of cavities into the microwell holder, 100 µl of the diluted sera and ready to use Controls are each added to the corresponding wells, the position A1 remains empty (blank). One Negative Control and two Standard Controls are carried along. The plate is incubated in a humid chamber at 37 °C for 30 min. The bottom of the cavities should not be in touch with materials that conduct temperature well (metal or wet paper).

The Controls that correspond to the determinations (IgG or IgM) are to be used.

A1	Blank
B1	Negative Control
C1	Standard Control
D1	Standard Control
E1, F1	patient serum 1, 2, etc.

#### **Attention!**

**The ELISA plate should not be placed in a cold incubation container which is heated to 37 °C during the incubation. The container must already be adapted to 37 °C prior to the incubation.**

### 10.5. Washing

Decant or aspirate all Microwells into a waste container with a disinfectant. Ensure complete removal of the liquid from the Microwells by tapping the inverted plate onto absorbent paper. Then wash all wells 4 times with 300 µl of prepared Washing Buffer. Be sure to remove residual washing solution by firmly tapping the inverted Microwells on absorbent paper after single washing steps.

**If a microplate washer is used, take care that the washer is adjusted to the used microplate type.**

### 10.6. Second incubation

Add 100 µl of RIDASCREEN® SeroG LD or SeroM LD (*Anti-human-IgG-* or *Anti-human-IgM-Conjugate*) into the corresponding wells (including A1). Incubate the plate for 30 min at 37 °C in a humid chamber.

### 10.7. Washing

Wash 4 times according to step 10.5.

### 10.8. Third incubation

Add 100 µl of RIDASCREEN® SeroSC into each well. Incubate the plate for 30 min at 37 °C in a humid chamber. Following the incubation, the reaction is stopped by adding 100 µl of RIDASCREEN® SeroStopp to each well. After careful mixing (soft tapping on the edge of the plate) the absorbance is measured in a microplate reader at 450 nm (reference wavelength  $\geq$  620 nm). Zero adjustment is done against the blank (position A1).

**Remark:**

**To remove moisture, wipe the bottom of the microplate before measuring.**

### Summary of the test procedure

1. Bring all reagents to room temperature
2. Dilute the RIDASCREEN® SeroWP
3. Dilute the serum samples
4. Pipet 100 µl of the Standard Control, the Negative Control or the diluted samples into the Microwells; 30 minutes incubation at 37 °C in a humid chamber
5. Discard the incubate and wash 4 times with 300 µl of Washing Buffer
6. Add 100 µl of RIDASCREEN® SeroG LD or SeroM LD; 30 minutes incubation at 37 °C in a humid chamber
7. Discard the incubate and wash 4 times with 300 µl of Washing Buffer
8. Add 100 µl of RIDASCREEN® SeroSC; 30 minutes incubation at 37 °C in a humid chamber
9. After addition of 100 µl RIDASCREEN® SeroStopp spectrophotometric determination at 450 nm (reference wavelength  $\geq$  620 nm)

## 11. Analysis

Analysis of the test can be done in three different ways:

1. using the standard curve enclosed
2. using the table of values (see data sheet enclosed)
3. mathematically by the 4-parameter calculation or the  $\alpha$ -method

**The blank (A1) must be subtracted from all OD values prior to the analysis.**

### 11.1. Quality control

For the quality control, Standard Control (double determination) and Negative Control must be carried along with each test procedure. The test was carried out correctly if the mean absorbance value of the Standard Control at 450/620 nm corresponds to the value range of the data sheet enclosed. If the difference of both single measurements is larger than 20 %, the test must be repeated. The Negative Control must show an absorbance value at 450/620 nm  $<$  0.3 and must be negative in its validation.

### 11.2. Analysis with the standard curve

If the analysis is to be carried out with the standard curve, first the daily variation must be corrected by a double determination of the Standard Control (see also pt. 11.1.). The correction factor F is calculated with the given target value of the Standard Control and the actually measured Standard Control value. The Lot-dependent target value is given on the data sheet enclosed.

$$F = \frac{\text{target value of the Standard Control}}{\text{mean absorbance value of the Standard Control}}$$

All OD values of the samples are multiplied by the factor F. With these corrected values, the corresponding Unit value is read from the standard curve.

### 11.3. Analysis with the table of values

The mean absorbance value of the Standard Control (see also pt. 11.1.) determines the column of the value range in the attached value table, which is decisive for the actual measurement. Within this column the measured absorbance value of the sample is classed to the suited absorbance range and the corresponding titer in Units/ml is read in the second column from the left.

	Units/ml	range of standard control	
		0.87 - 0.92	
-	< 10.0	< 0.33	
?	10.0 - 14.0	0.33 - 0.40	
+	14.1 - 35.0	0.41 - 0.72	
	35.1 - 80.0	0.73 - 1.14	
	80.1 - 150.0	1.15 - 1.50	
	150.1 - 300.0	1.51 - 1.84	
	300.1 - 1,000	1.85 - 2.34	
	> 1,000	> 2.34	

Fig. 1: Example for IgG determination  
(Excerpt from a Lot dependent data sheet)

The mean absorbance value in a measurement is 0.89 for example. In this case, the column with the value range from 0.87 to 0.92 is decisive for result determination. Then, a patient sample with an absorbance of 0.48 lies in a titer range from 14 to 35 Units/ml. (The foregoing values are an example and can differ from the actual values of the Lot-dependent data sheet.)

The valuation of the determined results - positive (+), negative (-), equivocal (?) - can be taken from the first column of the table of values on the data sheet enclosed.

### 11.4. Mathematical analysis

Needed values for a mathematical analysis by the 4-parameter calculation or the  $\alpha$ -method are given on the data sheet enclosed.

### 11.5. Test result

Tab. 1: Valuation of calculated Units

	IgG	IgM
negative	< 10 U/ml	< 10 U/ml
equivocal	10 - 14 U/ml	10 - 13 U/ml
positive	> 14 U/ml	> 13 U/ml

## 12. Remarks about the test procedure and interpretation

The RIDASCREEN® EBV VCA EIA detects antibodies against capsid antigen of EBV. IgM antibodies against this antigen can be used as a early reference to an acute infection. Absence of IgM with a positive IgG result indicates for a latent EBV infection. For conformation of the latent stage, detection of IgG antibodies against nuclear antigen of EBV should be done (RIDASCREEN® EBV EBNA, Art. No. K6621). Those antibodies are only detectable with the transition of the virus into its latent stage.

In comparison to IgM, IgG antibodies against VCA are later detectable in the course of first infection. Then, they rise significantly and fall down slightly after having reached the latent stage. IgG titer remain detectable for the whole life.

IgG antibodies against EBV can be found in 95% of the population. As the virus remains latently in a part of the B cells after got over an acute infection, the antibody-titer against EBV could be increased again after a stimulation of the immune system; this has also an effect on the formation of IgM antibodies. Therefore, the achieved results have always to be interpreted in comparison with the clinical view. To improve the diagnostic results, it is advisable to examine two successively taken sera of one patient.

## Appendix

### Literature

1. Burkhardt, F.: Mikrobiologische Diagnostik. Georg Thieme Verlag, Stuttgart (1992)
2. Gorgievski-Hrisoho, M., et al.: Serodiagnosis of infectious mononucleosis by using recombinant Epstein-Barr virus antigens and enzyme-linked immunosorbent assay technology. *J. Clin. Mikrobiol.* 28: 2305-2311 (1990)
3. Henle, W., et al.: Antibody responses to Epstein-Barr virus-determined nuclear antigen (EBNA)-1 and EBNA-2 in acute and chronic Epstein-Barr virus infection. *Proc. Natl. Acad. Sci. USA* 84: 570-574 (1987)
4. Hille, A., et al.: Expression of Epstein-Barr virus nuclear antigen 1, 2A and 2B in the baculovirus expression system: serological evaluation of human antibodies to these proteins. *J. Med. Virol.* 39: 233-241 (1993)
5. Linde A.: Diagnosis of Epstein-Barr virus-related diseases. *Scand. J. Infect. Dis. Suppl.* 100: 83-8 (1996)
6. Milman, G., et al.: Carboxyl-terminal domain of the Epstein-Barr virus nuclear antigen is highly immunogenic in man. *Proc. Natl. Acad. Sci. USA* 82: 6300-6304 (1985)
7. Modrow, S., Falke, D.: *Molekulare Virologie*. Spektrum Akademischer Verlag, Heidelberg (1997)
8. Vroman, B., et al.: Characterization of a major protein with a molecular weight of 160,000 associated with the viral capsid of Epstein-Barr virus. *J. Virol.* 53: 107-113 (1985)