RIDASCREEN[®] Mycoplasma

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Enzyme immunoassay for the detection of antibodies
(IgA, IgG, IgM) against Mycoplasma pneumoniae in serum

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—————

RIDA® und RIDASCREEN®
are registered trademarks of R-Biopharm AG
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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Contents

	page
1. Intended use	17
2. General	17
3. Test principle	18
4. Reagents provided	19
5. Materials required but not provided	21
6. Warnings and precautions for the users	21
7. Storage instructions	22
8. Indication of instability or deterioration of reagents	22
9. Specimen collection and storage	22
10. Test procedure	23
11. Analysis	26
12. Remarks about the test procedure and interpretation	28

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RIDASCREEN[®] Mycoplasma

Enzyme immunoassay for the detection of antibodies (IgA, IgG, IgM) against *Mycoplasma pneumoniae* in serum

1. Intended use

The RIDASCREEN[®] Mycoplasma test is an in vitro diagnosticum for specific detection of IgA, IgG and IgM antibodies against *Mycoplasma pneumoniae* in human serum.

2. General

Mycoplasma are bacteria without cell walls which are characterized by their pleomorphism. They have a size of 0.4 µm in their coccoid form. In contact with host cells or inert material such as glass or plastic, they take on a spindle-like form with filamentous structures. The most important representative of human-pathogenic *Mycoplasma* of the respiratory tract is *Mycoplasma pneumoniae*.

The incubation period of an infection with *M. pneumoniae* takes a relatively long time and lasts from 10 days to three weeks. After the clinical symptoms have subsided, *M. pneumoniae* can be excreted for weeks. However, the infectivity of the causative agent is very slight, so that close contact must exist for transmission to take place, such as given in families, day care centers, schools and military facilities. Illnesses occur mainly during childhood and early adulthood. In young children under the age of five, the infection often occurs almost symptomless. First contact with the causative agent does not produce immunity. Second infections are often more severe. Over-shooting reactions of the immunologically competent host are made responsible.

The clinical picture of a *M. pneumoniae* infection is that of a tracheobronchitis with lethargy, fever, headache and persistent, initially unproductive coughing. In approx. 10 % of all cases, atypical, interstitial pneumonia develops. This development can be very difficult but seldom takes a deadly turn. In addition, a multitude of extrapulmonary complications can be found like in other bacterial and viral infections and which are not *M. pneumoniae* specific. These are among others: Ery-

thema multiforme, Erythema nodosum, arthralgia, myalgia, hemolytical anaemia, meningitis, polyneuritis, myocarditis und pericarditis. A differential diagnostical clarificaton is, therefore, very important.

The RIDASCREEN® Mycoplasma test is, as an enzyme immunoassay, a simple and highly sensitive method for the detection of antibodies against *M. pneumoniae*. In contrast to other immunological procedures such as the CFT, a differentiation of the different immunoglobulin classes is possible with the EIA.

3. Test principle

This test is an enzyme immunoassay (EIA) for the detection of antibodies. On the surface of the microtiter wells, purified antigens of *Mycoplasma pneumoniae* are bound. Diluted serum samples and controls are pipetted into the wells and incubated at 37 °C in a humid chamber. Present antibodies bind to the immobilized antigen. Unbound material is removed in a washing step.

In a second step, a peroxidase-conjugated anti-human antibody (anti-IgA, -IgG or -IgM) is added. After incubation, unbound conjugate is removed by washing. Substrate (H_2O_2/TMB) is added to the wells and incubated at 37 °C in a humid chamber. The enzyme bound in the wells converts the colorless substrate to a blue color. Addition of stop solution converts the color from blue to yellow. The absorption is measured at 450 nm wavelength (reference wavelength ≥ 620 nm).

4. Reagents provided

4.1. RIDASCREEN® Mycoplasma IgA (Art. No. K 4311)

The reagents in one package are sufficient for 96 determinations.
Each test kit contains:

- 1 x 12 Microwell Strips (breakable) with 8 wells each in a frame;
coated with *M. pneumoniae* antigen;
in a resealable foil bag
- 1 x RIDASCREEN® SeroPP (110 ml);
Sample Diluent, ready to use;
contains 0.01% Thimerosal
- 1 x RIDASCREEN® SeroWP (100 ml, 10x conc.);
Washing Buffer, contains 0.2 % Bronidox-L
- 1 x Standard Control (2.5 ml);
diluted human serum, ready to use;
contains 0.01% Thimerosal
- 1 x Negative Control (1.2 ml);
diluted human serum, ready to use;
contains 0.01% Thimerosal
- 1 x RIDASCREEN® SeroA LD (12 ml);
Anti-human-IgA-Conjugate, ready to use;
peroxidase-conjugated antibody (goat)
- 1 x RIDASCREEN® SeroSC (12 ml);
H₂O₂/tetramethylbenzidine (TMB), ready to use
- 1 x RIDASCREEN® SeroStopp (12 ml);
1 M sulfuric acid
- 1 x Instructions for use

4.2. RIDASCREEN® Mycoplasma IgG (Art. No. K 4321)

The reagents in one package are sufficient for 96 determinations.
Each test kit contains:

- 1 x 12 Microwell Strips (breakable) with 8 wells each in a frame;
coated with *M. pneumoniae* antigen;
in a resealable foil bag
- 1 x RIDASCREEN® SeroPP (110 ml);
Sample Diluent, ready to use;
contains 0.01% Thimerosal
- 1 x RIDASCREEN® SeroWP (100 ml, 10x conc.);
Washing Buffer, contains 0.2 % Bronidox-L

- 1 x Standard Control (2.5 ml);
diluted human serum, ready to use;
contains 0.01% Thimerosal
- 1 x Negative Control (1.2 ml);
diluted human serum, ready to use;
contains 0.01% Thimerosal
- 1 x RIDASCREEN[®] SeroG HD(12 ml);
Anti-*human-IgG*-Conjugate, ready to use;
peroxidase-conjugated antibody (rabbit)
- 1 x RIDASCREEN[®] SeroSC (12 ml);
H₂O₂/tetramethylbenzidine (TMB), ready to use
- 1 x RIDASCREEN[®] SeroStopp (12 ml);
1 M sulfuric acid
- 1 x Instructions for use

4.3. RIDASCREEN[®] Mycoplasma IgM (Art. No. K 4331)

The reagents in one package are sufficient for 96 determinations.
Each test kit contains:

- 1 x 12 Microwell Strips (breakable) with 8 wells each in a frame;
coated with *M. pneumoniae* antigen;
in a resealable foil bag
- 1 x RIDASCREEN[®] SeroPP (110 ml);
Sample Diluent, ready to use;
contains 0.01% Thimerosal
- 1 x RIDASCREEN[®] SeroWP (100 ml, 10x conc.);
Washing Buffer, contains 0.2 % Bronidox-L
- 1 x Standard Control (2.5 ml);
diluted human serum, ready to use;
contains 0.01% Thimerosal
- 1 x Negative Control (1.2 ml);
diluted human serum, ready to use;
contains 0.01% Thimerosal
- 1 x RIDASCREEN[®] SeroM HD (12 ml);
Anti-*human-IgM*-Conjugate, ready to use;
peroxidase-conjugated antibody (rabbit)
- 1 x RIDASCREEN[®] SeroSC (12 ml);
H₂O₂/tetramethylbenzidine (TMB), ready to use
- 1 x RIDASCREEN[®] SeroStopp (12 ml);
1 M sulfuric acid
- 1 x Instructions for use

5. Reagents required but not provided

5.1. Reagents

- Distilled or deionized water

5.2. Accessories

- Moist chamber at 37 °C
- Test tubes
- Vortex mixer
- Micropipets for volumes of 10 - 100 µl and 100 - 1000 µl
- Measuring cylinder (1000 ml)
- Microplate washer or multichannel pipet
- Microplate reader (450 nm, reference wavelength \geq 620 nm)
- Absorbent paper

6. Warnings and precautions for the users

The control sera (standard control and negative control) have been tested for HIV- and HCV-Ab as well as for HBsAg and were found to be negative. However, they as well as the patient samples should be considered potentially contagious and be treated with the necessary safety precautions.

The standard control and the negative control as well as the sample diluent contain 0.01 % Thimerosal. Contact with skin or mucous membranes must be avoided.

The washing buffer contains Bronidox-L. Contact with skin or mucous membranes must be avoided.

Hydrogen peroxide can cause cauterization. Handle with care!

The stop solution contains 1 N sulfuric acid. Avoid contact with skin and clothing!

All reagents and materials coming in contact with potential infectious specimens must be treated with disinfectants or autoclaved at 121 °C for at least one hour.

An exchange of the controls between kits of different lot numbers is not possible.

Microwell strips and reagents must not be used if pouch is damaged or vials are leaking.

7. Storage instructions

All reagents have to be stored at 2 - 8 °C and can be used up to the expiry date printed on the labels. Microbial contamination has to be avoided. A quality warranty cannot be given beyond the kit expiration date.

The diluted washing buffer (RIDASCREEN® SeroWP) can be used up to the expiry date printed on the labels if stored at 2 - 8 °C.

Allow reagents and microwell strips to get room temperature before use. To avoid moisture within the strips, do not take the strips out of the foil bag before having reached room temperature. The foil bag should be opened with a pair of scissors without detaching the fastener. Return any unused strips to the foil bag, reseal and store them directly at 2 - 8 °C.

It is absolutely necessary to avoid contamination of the substrate (RIDASCREEN® SeroSC) with the conjugate solution, since this results in a coloration of the substrate. In the same way, the colorless substrate must be protected from exposure to direct light to avoid deterioration or coloration by autooxidation. If the substrate has turned blue, the reagent should be discarded.

8. Indication of instability or deterioration of reagents

The following criteria may indicate a reagent deterioration:

- a turbidity or a blue coloration of the colorless substrate prior to its use
- an absorbance value of the negative control at 450/620 nm > 0.3
- an absorbance value of the positive control at 450/620 nm out of the values given in the attached tables of values

9. Specimen collection and storage

The RIDASCREEN® Mycoplasma EIA has been evaluated for the investigation of human serum samples. Repeated freezing and thawing of the samples as well as microbial contamination must be avoided. The application of heat treated, lipemic, hemolytic, icteric or turbid samples can lead to wrong results.

The sample material can be stored for up to 1 week at 2 - 8 °C if the test cannot be carried out immediately. A prolonged storage of samples is possible at -20 °C.

Diluted samples can be used for up to 7 hours if they have been stored at 2 - 8 °C.

10. Test procedure

10.1. Preliminary comments

The test has to be used only by experienced laboratory personal. Please refer to guidelines for safety regulations in medical laboratories. The test protocol must be followed strictly.

The RIDASCREEN® EIAs for serological antibody detection contain some test independent reagents. These reagents are labeled RIDASCREEN® Sero (SeroWP, SeroPP, SeroA LD, SeroG HD, SeroM HD, SeroSC, SeroStopp). They are not kit specific and can be changed between single test kits with correspondingly labeled reagents. But, this is not valid for the controls.

Bring all reagents and the microwell strips to room temperature before use. Mix the reagents well before use. Reproducibility in any EIA depends on exact pipetting, the observance of incubation times and temperature and the consistency of wash sequences. A deviation from the given incubation times and temperatures leads to an excessive difference of the standard control serum to the given target value, which is no longer covered by the determined area of values. During the washing steps, take care that all wells are filled with buffer and that the liquid is completely removed from the wells. Do not allow microwells to dry between steps. Avoid direct sunlight during all incubations. Covering the microtiter plate is recommended.

Except the washing buffer, all reagents are ready to use.

10.2. Preparation of the washing buffer

1 part of the concentrated washing buffer (RIDASCREEN® SeroWP) is diluted with 9 parts of distilled water. Crystals in the buffer concentrate can be dissolved in a waterbath at 37 °C. Add 100 ml of the concentrated washing buffer to a 1000 ml graduated cylinder. Bring the final volume to 1000 ml with distilled or deionized water. The diluted washing buffer can be used up to the expiry date printed on the labels if stored at 2 - 8 °C.

10.3. Preparation of the samples

Before starting the test, serum samples have to be diluted 1:100 with the sample diluent (RIDASCREEN® SeroPP).

Dilution of serum samples (IgA, IgG and IgM determination):

e.g. 10 µl serum + 990 µl sample diluent

For IgM determination, sera should be treated with a IgG-absorbent (e.g. RIDA® RF-Absorbens, Art. No. Z 0202) prior to the test, following that it should be brought to its final dilution of 1:100 with the sample diluent.

Attention!

Negative control and standard control are ready to use and must not be diluted or absorbed.

10.4. First incubation

After insertion of a sufficient number of cavities into the microwell holder, 100 µl of the diluted sera and ready to use controls are each added to the corresponding wells, the position A1 remains empty (blank). One negative control and two standard controls are carried along. The plate is incubated in a humid chamber at 37 °C for 30 min. The bottom of the cavities should not be in touch with materials that conduct temperature well (metal or wet paper).

The controls that correspond to the determinations (IgA, IgG or IgM) are to be used.

A1	blank
B1	negative control
C1	standard control
D1	standard control
E1, F1	patient serum 1, 2, etc.

Attention!

The ELISA plate should not be placed in a cold incubation container which is heated to 37 °C during the incubation. The container must already be adapted to 37 °C prior to the incubation.

10.5. Washing

Decant or aspirate all microwells into a waste container with a disinfectant. Ensure complete removal of the liquid from the microwells by tapping the inverted plate onto absorbent paper. Then wash all wells 4 times with 300 µl of prepared washing buffer. Be sure to remove residual washing solution by firmly tapping the inverted microwells on absorbent paper after single washing steps.

If a microplate washer is used, take care that the washer is adjusted to the used microplate type.

10.6. Second incubation

Add 100 µl of RIDASCREEN® SeroA LD, SeroG HD or SeroM HD (anti-*human-IgA*-, anti-*human-IgG*- or anti-*human-IgM*-conjugate) into the corresponding wells (including A1). Incubate the plate for 30 min at 37 °C in a humid chamber.

10.7. Washing

Wash 4 times according to step 10.5.

10.8. Third incubation

Add 100 µl of RIDASCREEN® SeroSC into each well. Incubate the plate for 30 min at 37 °C in a humid chamber. Following the incubation, the reaction is stopped by adding 100 µl of RIDASCREEN® SeroStopp to each well. After careful mixing (soft tapping on the edge of the plate) the absorbance is measured in a microplate reader at 450 nm (reference wavelength \geq 620 nm). Zero adjustment is done against the blank (position A1).

Remark:

To remove moisture, wipe the bottom of the microplate before measuring.

Summary of the test procedure

1. Bring all reagents to room temperature
2. Dilute the RIDASCREEN® SeroWP
3. Dilute the serum samples
4. Pipet 100 µl of the standard control, the negative control or the diluted samples into the microwells; 30 minutes incubation at 37 °C in a humid chamber
5. Discard the incubate and wash 4 times with 300 µl of washing buffer
6. Add 100 µl of RIDASCREEN® SeroA LD, SeroG HD or SeroM HD; 30 minutes incubation at 37 °C in a humid chamber
7. Discard the incubate and wash 4 times with 300 µl of washing buffer
8. Add 100 µl of RIDASCREEN® SeroSC; 30 minutes incubation at 37 °C in a humid chamber
9. After addition of 100 µl RIDASCREEN® SeroStopp spectrophotometric determination at 450 nm (reference wavelength \geq 620 nm)

11. Analysis

Analysis of the test can be done in three different ways:

1. using the standard curve enclosed
2. using the table of values (see data sheet enclosed)
3. mathematically by the 4-parameter calculation or the α -method

The blank (A1) must be subtracted from all OD values prior to the analysis.

11.1. Quality control

For the quality control, standard control (double determination) and negative control must be carried along with each test procedure. The test was carried out correctly if the mean absorbance value of the standard control at 450/620 nm corresponds to the value range of the data sheet enclosed. If the single measurements diverge from the mean absorbance value more than 20 %, the test must be repeated. The negative control must show an absorbance value at 450/620 nm < 0.3.

If the expected control values are not fulfilled, please check the following before repeating the test:

- expiration date of the reagents
- calibration of the used instruments
- exact test procedure
- visual examination of kit components for signs of contamination, deterioration or leakage; substrate solution must not be used if turned blue

If the control data are not fulfilled after repeating, please contact your local distributor of R-Biopharm.

11.2. Analysis with the standard curve

If the analysis is to be carried out with the standard curve, first the daily variation must be corrected by a double determination of the standard control (see also pt. 11.1.). The correction factor F is calculated with the given nominal value of the standard control and the actually measured standard control value. The Lot-dependent nominal value is given on the data sheet enclosed.

$$F = \frac{\text{nominal value of the standard control}}{\text{mean absorbance value of the standard control}}$$

All OD values of the samples are multiplied by the factor F. With these corrected values, the corresponding Unit value is read from the standard curve.

11.3. Analysis with the table of values

The mean absorbance value of the standard control (see also pt. 11.1.) determines the column of the value range in the attached value table, which is decisive for the actual measurement. Within this column the measured absorbance value of the sample is classed to the suited absorbance range and the corresponding titer in Units/ml is read in the second column from the left.

The mean absorbance value of the standard control in a measurement is 0.86 for example. In this case, the column with the value range from 0.84 to 0.88 is decisive for result determination. Then, a patient sample with an absorbance of 0.65 lies in a titer range from 31.1 to 70.0 Units/ml. (The foregoing values are an example and can differ from the actual values of the Lot-dependent data sheet.)

The valuation of the determined results - positive (+), negative (-), equivocal (?) - can be taken from the first column of the table of values on the data sheet enclosed.

	U/ml	Range of standard control	
			0.84 - 0.88
-	< 23.0		< 0.30
?	23.0 - 31.0		0.30 - 0.38
+	31.1 - 70.0		0.39 - 0.69
	70.1 - 120.0		0.70 - 0.98
	120.1 - 200.0		0.99 - 1.32
	200.1 - 400.0		1.33 - 1.86
	400.1 - 1000.0		1.87 - 2.69
	> 1000.0		> 2.69

Fig. 1: Example for IgG determination
(Excerpt from a Lot dependent data sheet)

11.4. Mathematical analysis

Needed values for a mathematical analysis by the 4-parameter calculation or the α -method are given on the data sheet enclosed.

11.5. Test result

Tab. 1: Valuation of calculated units

	IgA	IgG	IgM
negativ	< 39 U/ml	< 23 U/ml	< 50 U/ml
grenzwertig	39 - 50 U/ml	23 - 31 U/ml	50 - 71 U/ml
positiv	> 50 U/ml	> 31 U/ml	> 71 U/ml

12. Remarks about the test procedure and interpretation

The RIDASCREEN[®] Mycoplasma EIA detects antibodies against *Mycoplasma pneumoniae*. It should be carried out if a well-founded suspicion of a *Mycoplasma* infection exists. A relation between the measured OD-value and the clinical relevance is not given. The results obtained should always be interpreted in connection with the clinical picture.

A negative result does not exclude a *Mycoplasma* infection, since the time at which the serum was taken could have been too early to detect any antibodies.

A positive IgG detection alone does not allow conclusions to be drawn about an acute illness. Important for the interpretation of an IgG diagnosis is the height and the development of the titer. Therefore, to improve the diagnostic statement, two consecutive sera of a patient should be examined. The results obtained should always be interpreted in connection with the clinical picture.

IgA antibodies allow conclusions to be drawn about a possible colonization of the mucous membrane with *M. pneumoniae*. Therefore, a high IgA titer indicates an acute infection. In some cases, IgA titers can also persist over a long period of time.

Usually, a positive IgM finding is only found after the first contact with the pathogen. But persistent IgM titers have also been observed.

Generally, two successive sera of a patient should always be examined to improve the diagnostic assertion.

A positive result does not exclude the presence of other pathogens as cause for an illness.

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